



### **Product Advantage: skin rejuvenation**

\*IPL handle use UK lamp, 150,000-30,000 shots. Spot size 15\*50mm. Cooling performance is excellent.

Ipl epilasyon makinesi tedarikçisi ipl gençleştirme makinesi ipl epilasyon ve cilt gençleştirme cihazı

\*Add your own logo and language into the machine.

### **Products Description**

# IPL+Laser

## Multifunctional Machine

IPL handle use UK lamp, 150,000-30,000 shots.  
Spot size 15\*50mm. Cooling performance is excellent.



420nm	510nm	560nm	585nm	640nm	690nm	755nm
Acne removal Vascular treatment	Pigmentation therapy, Remove freckles	Skin rejuvenation wrinkle removal	Remove Red blood silk	Hair removal (light skin)	Hair removal (brown skin)	Hair removal (black skin)

#### Color touch screenskin rejuvenation,One-button startup

1. 15.6 inch big touch screen, humanized designed
2. Male/femal, body parts, 6 skin types for option.
3. Englishyour native language.
4. Free OEM program service(logo/system...)

# SMART SYSTEM

Full screen, Wide field of view

15.6" HD SCREEN





EXFU

# One touch power on

ipl hair removal and skin rejuvenation device  
ipl hair removal and skin rejuvenation new  
ipl hair removal skin rejuvenation machine



ipl hair removal and skin rejuvenation machine

ipl hair removal skin rejuvenation machine

ipl rejuvenation machine

**24-hour online service!Welcome to consult!**